

ICEROSS SEAL-IN® X5 TRANSFEMORAL LINER




The Iceross Seal-In X5 transfemoral liner is the first liner of its kind with five integrated seals that deliver the most comfortable and secure suspension to date.


Its airtight fit is thanks to the flexible silicone suction suspension system, and its minimized pistoning and enhanced rotational control allows patients to enjoy greater comfort and stability.

It is recommended that Iceross Seal-In X5 Transfemoral Liner is used in conjunction with the Icelock® 500 Series Valve.

USER INFORMATION	
Amputation Level:	Transfemoral
Impact Level:	All
LINER INFORMATION	
Size Standard:	25,26.5,28,30,32,34,36,38,40,45,50,55.
Size Conical:	25,26.5,28,30,32,34,36,38.
Profile:	3mm
Matrix:	Full Length



STANDARD



CONICAL



ICEROSS®
Transfemoral

SUSPENSION METHOD



SEAL-IN X5

FEATURES



MATRIX



FABRIC
COVER



ACTIVE SKIN
CARE



MENTHOL



SILKEN INNER
SURFACE



CONICAL

PART NUMBER INFORMATION - ICEROSS SEAL-IN X5 TRANSFEMORAL LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-TF673XX	Standard 3mm	Seal-In X5	Full Length
I-TF678XX	Conical 3mm	Seal-In X5	Full Length